

# Ipas MVA Plus<sup>®</sup>

## For uterine evacuation and endometrial biopsy

Ipas MVA Plus<sup>®</sup> is a safe, effective and affordable medical device used for uterine evacuation up to 12 weeks gestation and endometrial biopsy. Manual Vacuum Aspiration (MVA) is the method of choice for uterine evacuation recommended by the World Health Organization (WHO) and the International Federation of Gynecology and Obstetrics (FIGO).



Product Code: PLUS

## PRODUCT DESCRIPTION



Reusable up to 25 times  
(where regulations allow)



Cylinder marks in 10 cc increments to  
measure evacuated uterine contents



Provides a minimum vacuum of 22  
inches or 558.8mm which is maintained  
for at least 30 minutes



Provides a holding capacity of 60cc



Sterilization<sup>1</sup> and High Level  
Disinfection<sup>2</sup> processing methods



Easy disassembly and reassembly



Latex-free plastic

## INSTRUCTION FOR USE

- Prepare the patient and provide them with pain relief - oral analgesia, paracervical block.
- Create a vacuum through the aspirator by pressing the valve buttons and pulling the plunger.
- Insert the cannula into the uterus and attach the prepared aspirator.
- Release the valve buttons on the aspirator and gently rotate to aspirate the uterine contents.

Scan Now: A  
Comprehensive Guide to  
the MVA Procedure



## INDICATIONS FOR USE

- Miscarriage management
- Induced or incomplete abortion up to 12 weeks since the woman's last menstrual period
- Endometrial biopsy

## COMPATIBILITY WITH CANNULAE

- Compatible with Ipas EasyGrip<sup>®</sup> Cannulae and size 12mm Flexible Karman Cannula
- Compatible with size 4-10mm Flexible Karman Cannulae with adapter
- Compatible with Endometrial Biopsy 3mm Cannula with a 6mm adapter

## BENEFITS

### Simple Solution

- Small, inexpensive, portable and quiet

### Easy access

- Ideal for performing procedures in the outpatient setting

### Convenient

- Provides the option of a procedure with local anaesthesia (paracervical block) with faster recovery time

### Suitable

- Ideal for settings with no electricity or with space constraints

Recommended Processing Methods for Ipas MVA Plus<sup>®</sup> - for further details, check the Instruction For Use (IFU) document

1. Sterilized by steam autoclavable for 30 minutes at 250°F (121°C); or by standard cold methods, disassembled must be fully immersed: Cidex<sup>®</sup> (2% glutaraldehyde solution) for 10 hours, Sporox<sup>®</sup> II for 6 hours.

2. High level disinfected in boiling water for 20 minutes; or by standard cold methods, disassembled parts must be fully immersed: Cidex<sup>®</sup> (2% glutaraldehyde solution) for 20 minutes, 0.5% Chlorine solution for 20 minutes, Sporox<sup>®</sup> II for 30 minutes.