

Tapered Uterine Dilators

For cervical dilation

Tapered Uterine Dilators are designed with tapered ends to be more gentle than non-tapered ends. Fully autoclavable for multiple uses.



Product Code: DD5-S

PRODUCT DESCRIPTION



Reusable up to 25 times (where regulations allow)



Sterilely packed



Set of 5 double-ended, tapered dilators: 5/6 mm, 7/8 mm, 9/10 mm, 11/12 mm, and 13/14 mm diameters



Approximately 29 cm (11.5 in) in length



Sterilization¹ and High Level Disinfection² processing methods



Latex-free, polypropylene plastic

INSTRUCTION FOR USE

- Introduce a dilator that has a diameter smaller than the estimated gestation of the pregnancy gently through the cervical canal, rotating the dilator with gentle pressure.
- Repeat this with progressively larger dilators until the internal cervical os has reached the appropriate for the patient's gestation.
- Remove dilator, insert cannula and start the aspiration procedure.

INDICATIONS FOR USE

For cervical dilation prior to surgical uterine evacuation

BENEFITS

Convenience

- Uterine depth markers

Lightweight

- Tapered - requires half the force of a non-tapered dilator

Safe

- Packed sterile

Recommended Processing Methods for Tapered Uterine Dilators - for further details, check the Instruction For Use (IFU) document

1. Sterilized by steam autoclavable for 30 minutes at 250°F (121°C); or by standard cold methods, items must be fully immersed: Cidex® (2% glutaraldehyde solution) for 10 hours, Sporox®II for 6 hours.

2. High level disinfected in boiling water for 20 minutes; or by standard cold methods, items must be fully immersed: Cidex® (2% glutaraldehyde solution) for 20 minutes, 0.5% Chlorine solution for 20 minutes, Sporox®II for 30 minutes.