

Ipas Miscarriage Management

A specialized kit designed to facilitate safe and effective MVA procedure



Product Code: MMKSU

PRODUCT DESCRIPTION

KIT COMPOSITION:



1 Ipas Double Valve Aspirator, the classic Ipas design



4 Ipas EasyGrip® Cannulae (6, 7, 8, and 9 mm diameters)



3 Tapered Uterine Dilators (5/6, 7/8, 9/10 mm diameters)

INSTRUCTION FOR USE

- Providers select the cannula size based on the woman's gestational age
- Prepare patient and provide pain relief – oral analgesia, paracervical block
- Creating the vacuum on the aspirator by pressing valve buttons and pulling plunger
- Insert cannula into the uterus and attach prepared aspirator
- Release valve buttons on the aspirator and gently rotate to aspirate uterine contents

INDICATIONS FOR USE

- Miscarriage management
- Endometrial biopsy
- Induced or incomplete abortion up to 12 weeks since the woman's last menstrual period

BENEFITS

Practical

- Ideal for emergency departments, ambulatory surgery centers and office-based settings

Ready to use

- Fully sterilized, pre-packaged, single-use kit

Convenient

- Portable, quiet and easily stored

Simple Use

- No electricity required and easy to disassemble and reassemble