

THE IPAS MANUAL VACUUM ASPIRATION TECHNOLOGY PRODUCT LINE CATALOGUE

A COST-EFFECTIVE SOLUTION FOR UTERINE EVACUATION¹

- BACKED BY 40 YEARS OF GLOBAL MARKET USE AND IPAS' LONGSTANDING QUALITY REPUTATION
- FDA-LISTED AND DESIGNED AND MANUFACTURED CONFORMING TO ISO 13485
- SMALL, INEXPENSIVE, PORTABLE, QUIET, AND SAFE^{2,3,4}
- IDEAL FOR PERFORMING PROCEDURES IN OUTPATIENT SETTING
- IDEAL FOR SETTINGS WITH NO ELECTRICITY OR WITH SPACE CONSTRAINTS
- LIMITED INVESTMENT REQUIRED
- EASY TO USE

CLINICAL BENEFITS IN MANUAL UTERINE EVACUATION

	IPAS ASPIRATORS	EVA	SHARP CURETTAGE
SAFE	X	X	
LOWER INCIDENCE OF SIDE EFFECTS/COMPLICATIONS ⁵	X	X	
COST EFFECTIVE	X		X*
PORTABLE	X		X
QUIET	X		X
SINGLE USE	X		
REUSABLE	X	X	X
OUTPATIENT/NO ANESTHESIA ⁶	X	X	
QUICK RECOVERY TIME	X	X	
NO ELECTRICITY REQUIRED	X		X
REQUIRES LIMITED SPACE	X		X
REQUIRES LIMITED INVESTMENT	X		X

*When used without anesthesia

References:

1. Veruyil DAA, Crowther CA. (1993). Suction v. conventional curettage in complete abortion: a randomized clinical trial. S Afr Med J.83:13-15.
2. Blumenthal, Paul D and R.E. Remsburg. 1994. A time and cost analysis of the management of incomplete abortion with manual vacuum aspiration. Int J Gynecol Obstet, 1994, 45: 261-267.
3. De Jonge ETM, Pattinson, RC, Makin JD et al. (1994). Is ward evacuation for uncomplicated incomplete abortion under systematic analgesia safe and effective: a randomized clinical trial. S Afr Med J, 84:481-483.
4. Fawcus S, McIntyre J, Jewkes RK et al. (1997). Management of incomplete abortions at South African public hospitals. S Afr Med J, 87:438-442.
5. World Health Organization, Department of Reproductive Health and Research. Safe Abortion: Technical and policy guidance for health systems. Second edition. 2012. Annex 5, Recommendation 1, page 113.
6. Hale RW, Kobara TY, Sharma SD, Tsuei JJ, Gramlich EP, Nakayama RT. Office termination

WomanCare Global Trading CIC is a UK Registered Community Interest Company (No.07231920).
E-mail: contact@dktwomancare.org
www.dktwomancare.org

N-CAT-NPS-501-EN Rev. 00

PRODUCT PORTFOLIO

ASPIRATORS

- Ipas MVA Plus® Aspirator
- Ipas Double Valve Aspirator
- Ipas Single Valve Aspirator

CANNULAE

- Ipas EasyGrip® Cannulae
- Flexible Karman Cannulae

DILATORS AND ACCESSORIES

- Adapter sets, dilators and aspirator accessories also available

KITS

- Miscarriage Management Kit

TRAINING TOOLS

- Pelvic Model

WomanCare GLOBAL A COMPANY OF **dkt** INTERNATIONAL

An ISO 9001-certified company



BENEFITS

Discreet and Patient-Friendly

- Quiet procedure meets patient's need for privacy
- Outpatient procedure enables faster waiting/recovery times

Safe and Effective

- 98% effective with low complication rates
- A World Health Organization (WHO)-recommended method for safe uterine evacuation

Affordable and Versatile

- Reduces costs, material and human resources related to general anesthesia/sedation and/or operating room
- Low initial investment and per-use cost

Gentle suction

- Reduces risk of uterine perforation compared to unsafe methods like sharp curettage
- Preserves products of conception (POC) better than electric vacuum aspiration (EVA) - helps identify end of pregnancy and diagnose miscarriage causes

Ipas MVA Plus® Aspirator

For uterine evacuation

REUSABLE · DISCREET · SAFE · COST-EFFECTIVE



CLINICAL INDICATIONS

- Induced or incomplete abortion up to 12 weeks since the woman's last menstrual period
- Miscarriage management
- Endometrial biopsy

PRODUCT CHARACTERISTICS

- Reusable up to 25 times (where regulations allow)
- Provides holding capacity of 60 cc
- Provides a minimum vacuum of 22 inches
- Maintains vacuum for at least 30 minutes
- Cylinder marks in 10 cc increments to measure evacuated uterine contents
- Compatible with Ipas EasyGrip® Cannulae and the 12 mm Flexible Karman Cannula
- Sterilized by steam autoclavable at 250°F (121°C); or by standard cold methods (Cidex®, Sporox®II)
- High level disinfected in boiling water; or by standard cold methods (Cidex®, 0.5% Chlorine, Sporox®II)
- Latex-free plastic
- Shelf life: 5 years from month of manufacture

ALSO AVAILABLE FROM WOMANCARE GLOBAL




IPAS DOUBLE VALVE ASPIRATOR

For uterine evacuation
Product Code: DVS-10 or DVS-SU



DENNISTON DILATORS

For cervical dilation associated with uterine evacuation
Product Code: DD5-S



IPAS SINGLE VALVE ASPIRATOR

For uterine evacuation
Product Code: SVS-10 or SVS-LF

Ipas Double-Valve Aspirator

For uterine evacuation and endometrial biopsy

EFFECTIVE · DISPOSABLE · AFFORDABLE · VERSATILE



CLINICAL INDICATIONS

- Induced or incomplete abortion up to 12 weeks since the woman's last menstrual period
- Miscarriage management
- Endometrial biopsy

PRODUCT CHARACTERISTICS

- 10 units sterilely packed
- Individually wrapped
- Ready to use
- Provides a holding capacity of 60 cc
- Provides a minimum vacuum of 24 inches
- Maintains vacuum for a minimum of 30 minutes
- Cylinder marks in 10 cc increments to measure evacuated uterine contents
- Compatible with Ipas EasyGrip® Cannulae and the 12 mm Flexible Karman Cannula
- Latex-free plastic
- Shelf life: 3 years from month of manufacture


ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS SINGLE VALVE ASPIRATOR

For uterine evacuation

Product Code: SVS-10 or SVS-LF



3MM CANNULA

For endometrial biopsy

Product Code: EB3



FLEXIBLE KARMAN CANNULAE

For uterine evacuation

Product Code: FC4 to FC10, and FC12

Ipas Double-Valve Aspirator

For uterine evacuation and endometrial biopsy

EFFECTIVE · DISPOSABLE · AFFORDABLE · VERSATILE



CLINICAL INDICATIONS

- Induced or incomplete abortion up to 12 weeks since the woman's last menstrual period
- Miscarriage management
- Endometrial biopsy

PRODUCT CHARACTERISTICS

- Individually packed
- Provides a holding capacity of 60 cc
- Provides a minimum vacuum of 24 inches
- Maintains vacuum for a minimum of 30 minutes
- Cylinder marks in 10 cc increments to measure evacuated uterine contents
- Compatible with Ipas EasyGrip® Cannulae and the 12 mm Flexible Karman Cannula
- Latex-free plastic
- Shelf life: 5 years from month of manufacture

ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS MVA PLUS® ASPIRATOR

For uterine evacuation
Product Code: PLUS



IPAS EASYGRIP® CANNULAE

For uterine evacuation
Product Code: SR4 to SR10, and SR12



FLEXIBLE KARMAN CANNULAE

For uterine evacuation
Product Code: FC4 to FC10, and FC12

ipas Single Valve Aspirator

For uterine evacuation and endometrial biopsy

EFFECTIVE · REUSABLE · AFFORDABLE · VERSATILE



CLINICAL INDICATIONS

- Induced or incomplete abortion up to 6 weeks since the woman's last menstrual period
- Miscarriage management
- Endometrial biopsy

PRODUCT CHARACTERISTICS


- 10 units sterilely packed
- Reusable up to 25 times (where regulations allow)
- Provides a holding capacity of 60 cc
- Provides a minimum vacuum of 24 inches
- Maintains vacuum for a minimum of 30 minutes
- Cylinder marks in 10 cc increments to measure evacuated uterine contents
- Compatible with Flexible Karman Cannula sizes 4-6 mm and the Endometrial Biopsy 3MM cannula
- Latex-free plastic
- Shelf life: 3 years from month of manufacture

ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS DOUBLE VALVE ASPIRATOR

For uterine evacuation
Product Code: DVS-10 or DVS-SU



3MM CANNULA

For endometrial biopsy
Product Code: EB3



DENNISTON DILATORS

For cervical dilation associated with uterine evacuation
Product Code: DD5-S

IPAS Single Valve Aspirator

For uterine evacuation and endometrial biopsy

EFFECTIVE · REUSABLE · AFFORDABLE · VERSATILE



CLINICAL INDICATIONS

- Induced or incomplete abortion up to 6 weeks since the woman's last menstrual period
- Miscarriage management
- Endometrial biopsy

PRODUCT CHARACTERISTICS

- Individually packed
- Reusable up to 25 times (where regulations allow)
- Provides a holding capacity of 60 cc
- Provides a minimum vacuum of 24 inches
- Maintains vacuum for a minimum of 30 minutes
- Cylinder marks in 10 cc increments to measure evacuated uterine contents
- Compatible with Flexible Karman Cannula sizes 4-6 mm and the Endometrial Biopsy 3MM cannula
- Latex-free plastic
- Shelf life: 3 years from month of manufacture

ALSO AVAILABLE FROM WOMANCARE GLOBAL



UTERINE PELVIC MODEL

For uterine evacuation training
Product Code: UPM



IPAS EASYGRIP® CANNULAE

For uterine evacuation
Product Code: SR4 to SR10, and SR12



FLEXIBLE KARMAN CANNULAE

For uterine evacuation
Product Code: FC4 to FC10, and FC12

IPAS EasyGrip® Cannulae

The tactile response of a rigid curette with the gentle probe of a flexible cannulae

SEMI-RIGID · INNOVATIVE · VERSATILE · REUSABLE



PRODUCT CHARACTERISTICS

- Ready to use and sterile packed
- Reusable up to 25 times
- Permanently-integrated winged base
- Color-coded by size
- Uterine depth markers from 6 to 11 cm printed on cannulae
- Sized 4 mm to 8 mm - two "whistle cut" openings
- Sized 9 mm, 10 mm, and 12 mm - one large "scoop cut" opening
- Approximately 24 cm (9.5 in) in length
- Latex-free plastic
- Shelf life: 5 years from month of manufacture

COLOR CODE



PRODUCT CODE/DIAMETER CHART

Size	4 mm	5 mm	6 mm	7 mm	8 mm	9 mm	10 mm	12 mm
Code	SR4	SR5	SR6	SR7	SR8	SR9	SR10	SR12

ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS DOUBLE VALVE ASPIRATOR

For uterine evacuation
Product Code: DVS-10 or DVS-SU



IPAS MVA PLUS® ASPIRATOR

For uterine evacuation
Product Code: PLUS



UTERINE PELVIC MODEL

For uterine evacuation training
Product Code: UPM

Flexible Karman Cannulae

The original flexible cannulae

EFFECTIVE · DISPOSABLE · AFFORDABLE · FLEXIBLE



PRODUCT CHARACTERISTICS

- Sterilely packed
- Printing on cannulae color-coded by size
- Uterine depth markers from 6 to 11 cm printed on cannulae
- Sized 4 mm to 8 mm - two "whistle cut" openings
- Sized 9 mm, 10 mm, and 12 mm - one large "scoop cut" opening
- Approximately 24 cm (9.5 in) in length
- Sizes 4 mm - 6 mm compatible with Ipas Single Valve Aspirator
- Adapters required for use with Ipas MVA Plus and Double Valve Aspirator
- Latex-free plastic
- Shelf life: 5 years from month of manufacture

COLOR CODE



PRODUCT CODE/DIAMETER CHART

Size	4 mm	5 mm	6 mm	7 mm	8 mm	9 mm	10 mm	12 mm
Code	FC4	FC5	FC6	FC7	FC8	FC9	FC10	FC12

WomanCare Global Trading CIC is a UK Registered Community Interest Company (No.07231920).
 Email: contact@dktwomancare.org
www.dktwomancare.org

©WomanCare Global
 N-PRO-FC#*-701-EN Rev. 01

ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS MVA PLUS® ASPIRATOR
 For uterine evacuation
 Product Code: PLUS



DENNISTON DILATORS
 For cervical dilation associated with uterine evacuation
 Product Code: DD5-S




IPAS SINGLE VALVE ASPIRATOR
 For uterine evacuation
 Product Code: SVS-10 or SVS-LF

3MM Cannula

For endometrial biopsy

EFFECTIVE · DISPOSABLE · AFFORDABLE · VERSATILE



CLINICAL APPLICATIONS

- Evaluation of infertility
- Evaluation of amenorrhea
- Evaluation of suspected endometrial cancer
- Evaluation of abnormal uterine bleeding
- Diagnosis of endometrial infections

PRODUCT CHARACTERISTICS

- Disposable and sterile packed
- Permanently integrated winged base
- Two opposing apertures
- Approximately 23 cm (9 in) in length
- Depth markings printed in cm
- Compatible with Ipas MVA Plus® Aspirator with a 6 mm adaptor
- Compatible with Ipas MVA Single Valve Aspirator
- Medical grade, latex-free copolyester plastic
- Shelf life: 5 years from month of manufacture

COLOR CODE

3mm

WomanCare Global Trading CIC is a UK Registered Community Interest Company (No.07231920).
Email: contact@dktwomancare.org
www.dktwomancare.org

©WomanCare Global
N-PRO-EB3-701-EN Rev. 01

ALSO AVAILABLE FROM WOMANCARE GLOBAL



UTERINE PELVIC MODEL

For uterine evacuation training
Product Code: UPM



FLEXIBLE KARMAN CANNULAE

For uterine evacuation
Product Code: FC4 to FC10, and FC12



IPAS MVA PLUS® ASPIRATOR

For uterine evacuation
Product Code: PLUS

Denniston Dilators

For cervical dilation

REUSABLE · GRADUATED · LIGHTWEIGHT



CLINICAL INDICATION

- For cervical dilation prior to uterine evacuation

PRODUCT CHARACTERISTICS

- Packed sterile
- Tapered - requires half the force of a non-tapered dilator
- Set of 5 double-ended, tapered dilators: 5/6 mm, 7/8 mm, 9/10 mm, 11/12 mm, and 13/14 mm diameters
- Approximately 29 cm (11.5 in) in length
- Latex-free, polypropylene plastic
- Sterilized by steam autoclavable at 250°F (121°C) or by standard cold methods (Cidex®, Sporox®II)
- High level disinfected in boiling water; or by standard cold methods (Cidex®, 0.5% Chlorine, Sporox®II)
- Shelf life: 5 years from month of manufacture

ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS MVA PLUS® ASPIRATOR
For uterine evacuation
Product Code: PLUS



IPAS EASYGRIP® CANNULAE
For uterine evacuation
Product Code: SR4 to SR10, and SR12
or SR4-SU to SR10-SU, and SR12-SU



IPAS SINGLE VALVE ASPIRATOR
For uterine evacuation
Product Code: SVS-10 or SVS-LF

Uterine Pelvic Model

For effective training

SIMPLE TO USE · OPEN DESIGN · REALISTIC



PRODUCT CHARACTERISTICS

- Open design for ease of intrauterine procedure training, including uterine evacuation
- Simulation of both smooth and gritty tactile sensations during uterine evacuation
- Plexiglass base with suction cups to ensure stability
- Lightweight and portable
- Measures 25 x 25 x 18 cm (10 x 10 x 10 in) approximately



ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS DOUBLE VALVE ASPIRATOR

For uterine evacuation
Product Code: DVS-10 or DVS-SU



DENNISTON DILATORS

For cervical dilation associated with uterine evacuation
Product Code: DD5-S




IPAS SINGLE VALVE ASPIRATOR

For uterine evacuation
Product Code: SVS-10 or SVS-LF

Miscarriage Management

For early pregnancy loss management

COST-EFFECTIVE · CONVENIENT · EFFICIENT




PRODUCT CHARACTERISTICS

- Fully sterilized, pre-packaged, single-use kit
- Hand-held, simple, and portable
- Affordable alternative to surgical operating room procedures
- Ideal for emergency departments, ambulatory surgery centers and office-based setting

KIT COMPOSITION

- 1 Ipas Double Valve Aspirator, the classic Ipas design
- 4 Ipas EasyGrip® Cannulae (6, 7, 8, and 9 mm diameters)
- 3 Denniston Dilators (5/6, 7/8, 9/10 mm diameters)

ALSO AVAILABLE FROM WOMANCARE GLOBAL



3MM CANNULA

For endometrial biopsy
Product Code: EB3



IPAS EASYGRIP® CANNULAE

For uterine evacuation
Product Code: SR4 to SR10, and SR12



DENNISTON DILATORS

For cervical dilation associated with uterine evacuation
Product Code: DD5-S



WomanCare | A COMPANY OF

GLOBAL | **dkt** | INTERNATIONAL

**CREATING A WORLD WHERE WOMEN ARE EMPOWERED
TO LIVE THEIR LIVES TO THE FULLEST**

©Copyright WomanCare Global Trading CIC | contact@dktwomancare.org | www.dktwomancare.org

WomanCare Global
1640 Roanoke Blvd,
Salem, VA 24153, USA

WomanCare Global Trading CIC
One Bartholomew Close,
EC1A 7BL, London, UK

DKT WomanCare Global Services SAS
28 Boulevard Haussmann
75009 Paris, France