THE IPAS MANUAL VACUUM **ASPIRATION TECHNOLOGY** PRODUCT LINE CATALOGUE

A COST-EFFECTIVE SOLUTION FOR UTERINE EVACUATION

CLINICAL BENEFITS IN MANUAL UTERINE EVACUATION

	IPAS Aspirators	EVA	SHARP CURETTAGE
SAFE	Х	Х	
LOWER INCIDENCE OF Side effects/complications ⁵	х	Х	
COST EFFECTIVE	Х		Х*
PORTABLE	X		Х
QUIET	χ		Х
SINGLE USE	Х		
REUSABLE	Х	Х	Х
OUTPATIENT/NO ANESTHESIA ⁶	Х	X	
QUICK RECOVERY TIME	X	Х	
NO ELECTRICITY REQUIRED	Х		х
REQUIRES LIMITED SPACE	Х		х
REQUIRES LIMITED INVESTMENT	х		х

*When used without anesthesia

- References:

 1. Veruyl DAA, Crowther CA. (1993). Suction v. conventional curettage in complete abortion: a randomized clinical trial.

 S Afr Med J,83:13-15.

 2. Blumenthal, Paul D and R.E. Remsburg. 1994. A time and cost analysis of the management of incomplete abortion with manual vacuum aspiration. Int J Gynecol Obstet, 1994, 45: 261-267.

 3. De Jonge ETM, Pattinson, RC, Makin JD et al. (1994). Is ward evacuation for uncomplicated incomplete abortion under systematic analgesia safe and effective: a randomized clinical trial. S Afr Med J, 84:481-483.

 4. Fawcus S, McIntye J, Jewkes RK et al. (1997). Management of incomplete abortions at South African public hospitals. S Afr Med J, 87:438-442.

 5. World Health Organization, Department of Reproductive Health and Research. Safe Abortion: Technical and policy guidance for health systems. Second edition. 2012. Annex S, Recommendation 1, page 113.

 6. Hale RW, Kobara TY, Sharma SD, Tsuei JJ, Gramlich EP, Nakayama RT. Office termination

- BACKED BY 40 YEARS OF GLOBAL MARKET USE AND IPAS' **LONGSTANDING QUALITY REPUTATION**
- FDA-LISTED AND DESIGNED AND MANUFACTURED CONFORMING TO ISO 13485
- SMALL, INEXPENSIVE, PORTABLE, QUIET. AND SAFE^{2,3,4}
- IDEAL FOR PERFORMING PROCEDURES IN OUTPATIENT SETTING
- IDEAL FOR SETTINGS WITH NO **ELECTRICITY OR WITH SPACE CONSTRAINTS**
- LIMITED INVESTMENT REQUIRED
- EASY TO USE

PRODUCT PORTFOLIO

ASPIRATORS

- Ipas MVA Plus® Aspirator
- Ipas Double Valve Aspirator
- Ipas Single Valve Aspirator

CANNULAE

- Ipas EasyGrip® Cannulae
- Flexible Karman Cannulae

DILATORS AND ACCESSORIES

• Adapter sets, dilators and aspirator accessories also available

KITS

Miscarriage Management Kit

TRAINING TOOLS

• Pelvic Model



Uterine Pelvic Model

For effective training

SIMPLE TO USE · OPEN DESIGN · REALISTIC



PRODUCT CHARACTERISTICS

- Open design for ease of intrauterine procedure training, including uterine evacuation
- Simulation of both smooth and gritty tactile sensations during uterine evacuation
- Plexiglass base with suction cups to ensure stability
- Lightweight and portable
- Measures 25 x 25 x 18 cm (10 x 10 x 10 in) approximately



ALSO AVAILABLE FROM WOMANCARE GLOBAL



IPAS DOUBLE VALVE ASPIRATOR
For uterine evacuation

Product Code: DVS-10 or DVS-SU



DENNISTON DILATORS

For cervical dilation associated with uterine evacuation Product Code: DD5-S



IPAS SINGLE VALVE ASPIRATOR

For uterine evacuation Product Code: SVS-10 or SVS-LF

